



Friesian Horse Society, Inc. (FHS)

Instructions for Microchip Implantation

Please provide the following instructions to your veterinarian:

Scan the microchip prior to implantation and make sure the number on the stickers corresponds to the microchip in the package that will be implanted. This will also insure that the microchip is in good working order before implantation in the horse.

The microchip is implanted in the nuchal ligament on the left side of the neck perpendicular to the point of the shoulder.
After implantation the chip should be scanned and verified. Please make sure to record the number and place the sticker on the Microchip form.

If there are any questions please contact the FHS office!

FHS - Friesian Horse Society Inc.
8090 Granite Falls Dr.
Redmond, OR 97756

Phone: 402-298-734
FHS Text: 402-804-3409

E-mail: info@friesianhorsesociety.com